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


For the management of fractures of the tibial head and proximal shaft, available methods employ plates, double plates, fixators and intramedullary procedures and combinations of these. Depending on the procedure, there is high soft tissue stress, high bone stress due to denuding, risk of implant failure and axial instability with the risk of secondary dislocation. This can lead to the development of a pseudarthrosis, implant loosening and fracture and delayed wound healing with subsequent infection, followed by the need for a change of procedure. The angle-stable, anatomically shaped lateral tibial head plate (4.5 mm) is now available for a soft tissue-sparing internal fixation that is as minimally invasive as possible and nevertheless free from loss of reduction. The implant presented here enables a primary definitive procedure to be carried out with sparing of soft tissues and preservation of bone perfusion. A unilateral procedure is possible because of the angular stability. The screw direction is anatomically optimised. Conventional 4.5 mm standard instruments are compatible.

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- extra-articular tibial head fractures (AO classification 41 - A2 and A3)
- selected articular tibial head fractures (AO 41 - B1 and B3)
- comminuted fractures of the proximal tibia (AO 41 - C1-C3)

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Anatomical shaped right and left Setting, bending, curving

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Anatomically adapted optimised screw direction

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Stability of plate body with angle-stable screws

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Supine position with freely movable draped leg

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Lateral access to the tibial head sparing the periosteum and preserving the muscle attachments

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Reduction and treatment of the concomitant tibial head fracture (B and C injury) and/or advancement of the plate in the required length into the fibialis anterior compartment – image intensifier control

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Fixation of the plate to the joint block (previously reconstructed if necessary)

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Alignment of the tibial shaft with regard to rotation, varus-valgus relative to the joint block - image intensifier control

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Fixation of the plate to the aligned tibial shaft - reduction of the tibial shaft to the anatomical plate possible through tension screw

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Separate incision over the distal plate holes and placement of unidirectional angle-stable shaft screws using drill guide

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Image intensifier control in both planes

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Closure of fascial incision, meniscus refixation if necessary, drain, skin suture

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
- Splinted position until removal of drain (2nd day)
- Physiotherapy after removal of drain
- Range of movement or orthosis after associated intra-articular injury
- Discharge with floor contact permitted depending on fracture type and associated injury at least 6-8 weeks

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Cancellous screws, D 6,0 mm (conical head thread, self-tapping)

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© Copyright Königsee Implantate GmbH
OSTEOSYNTHESE GmbH
KÖNIGSEE IMPLANTATE UND INSTRUMENTE ZUR OSTEOSYNTHESE GmbH