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

Forefoot screws

- Titan
titanium

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- **Drill-Fix forefoot screw**
- **Fracture compressing screw**

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Joint-preserving foot operations with corrective osteotomies of the metatarsal bones have the aim of eliminating pain and preserving or restoring foot function which is as physiological as possible. They are usually performed when the intermetatarsal angle is enlarged ($>10^\circ$) and are combined with distal soft tissue procedures, especially lateral release on the first digital ray (tenotomy of the abductor hallucis muscle). However, additional release on the lateral side of the metatarsal head (especially in the region of the shortened lateral joint capsule) worsens the blood supply of the osteotomized head and increases the risk of underperfusion and thus of bone necrosis.

Measurement of the intermetatarsal angle can be limited because radiographs are frequently not exactly focussed and the films are taken without loading. Good orientation is offered by the position of the sesamoid bones relative to the head of the metatarsal. If the lateral sesamoid is displaced into the intermetatarsal space, it can be assumed that the intermetatarsal angle has been widened. The aim of corrective osteotomy is to position the metatarsal head over the sesamoid bones again. To achieve this a lateral shift of the medially displaced head is necessary. If the joint surface is intact, mobility is preserved in the metatarsophalangeal joint of the great toe, which permits rolling off the foot over the first ray after such operations and powerful pushing off the weightbearing bone through the great toe.

Excessively long and non-contracted claw toes II-V can be corrected by shortening osteotomies of the metatarsal bones.

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Soft tissue correction on the lateral side (lateral release)

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Contracted valgus deformity of the big toe (precedes the corrective osteotomy)

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Skin incision between the first and second ray for approx. 3 cm, commencing just proximal to the interdigital web fold, blunt dissection down to the tendon of abductor hallucis and its attachment is divided sharply at the base of the lateral proximal phalanx and the lateral sesamoid. In addition, the lateral capsule of the MP joint, which is often shortened by the deformity of the great toe, can be incised parallel to the joint line. However, this involves the risk of further interference with the blood supply of the head of the metatarsal and is therefore controversial. This is followed by manual stretching by pressure on the great toe in the varus direction.

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Hallux valgus with enlargement of the intermetatarsal angle, adequate corrected position with maximum displacement of the metatarsal head by one third of the width of the shaft.

Contraindication: Significant osteoarthritis in the MP joint of the great toe.

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
Skin incision on the medial side of the 1st digital ray, slightly more towards the extensor aspect, extending from the proximal phalanx over the metatarsal head. If necessary, bursectomy, longitudinal opening of the capsule (variant: opening the capsule by making a proximally based capsule flap approx. 1 cm wide). Exposure of the pseudo-exostosis and removal of this preserving the plantar sulcus in order to avoid the unwanted development of hallux varus. Marking of the centre of the head with a Kirschner wire, which is inserted from the medial side. Using this wire as a guide, V-shaped osteotomy, open proximally. The angle enclosed by the V should not be greater than 60 degrees. If the plantar part of the incision is longer, fixation is easier and this also allows the tip of the screw to be situated outside the plantar joint surface. Depending on the position of the V-shaped osteotomy, the length of the ray can also be altered and plantar/dorsal displacement of the metatarsal head can take place. For instance, an osteotomy in a proximal and lateral direction leads to shortening of the metatarsal. Displacement of the head by a maximum of one third of the width of the shaft with the aim of positioning the head over the sesamoid bones. Temporary fixation with a threaded wire through which the length of screw required is measured with the length gauge. Drill with the two-step drill as far as the opposite cortex on the plantar side and on the dorsum cut the countersink hole for the head. There are two stepped drills on the screen for this purpose (normally, the shorter one is sufficient). Then introduction of the appropriately measured fracture compression screw over the threaded wire which is still in place. The screw should exert compression on the osteotomy, the head should be countersunk and the opposite cortex should be gripped. Now remove the projecting bone edge creating a harmonious transition to the diaphysis of the metatarsal. Capsule closure doubling the capsule to correct the position of the great toe and drawing the sesamoids medially (variant: distal shortening of the proximally based capsule flap according to the correct position of the great toe with refixation distally, followed by capsule closure). Skin suture. Dressing to splint the forefoot to support the corrected position of the great toe.

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Mobilisation with forefoot-relieving sandal for 6 weeks. If necessary, hallux valgus splint at night or adhesive strapping of the great toe for 6 weeks. Check radiograph before progressing to full loading of the forefoot. The screw can be left in situ.

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Hallux valgus with enlargement of the intermetatarsal angle, adequate corrected position with maximum displacement of the metatarsal head by one half of the width of the shaft.

Contraindication: Significant osteoarthritis in the MP joint of the great toe.

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
Skin incision on the medial side of the 1st metatarsal, slightly more towards the extensor aspect, if necessary, bursectomy, longitudinal opening of the capsule, exposure of the pseudo-exostosis and removal of this preserving the plantar sulcus in order to avoid the unwanted development of hallux varus. Subperiosteal exposure of the metatarsals sparing the vascular soft tissues. Longitudinal incision with the oscillating saw through the medial cortex as a subtotal osteotomy, more plantar proximally and more dorsal distally. The two horizontal osteotomies are each made at an angle of 45 degrees to the longitudinal incision, plantar proximally and dorsal distally. Then completion of the longitudinal osteotomy dividing the lateral cortex in the same plane as where it commenced medially. Mobilisation and lateral displacement of the plantar part with the undamaged metatarsal head by up to two thirds of the width of the shaft. Temporary fixation with a bone-holding forceps, then introduction of two threaded wires vertically from the dorsal aspect as far as the opposite plantar cortex (variant: the distal threaded wire can optionally be directed obliquely in the direction of the metatarsal head and should then lead the opposite cortex untouched). Measurement of the length with the length gauge over the introduced threaded wires. Drill with the two-step drill as far as the opposite cortex on the plantar side (without touching the opposite cortex in the variant with oblique insertion of the wire) and on the dorsum cut the countersink hole for the head. There are two stepped drills on the screen for this purpose. Then introduction of the two appropriately measured fracture compression screws over the threaded wires which are still in place. The screws should exert compression on the osteotomy, the head should be countersunk and the opposite cortex should be gripped (in the variant with the screw inserted obliquely towards the metatarsal head, leave the opposite cortex untouched). Now remove the distally projecting bone edge to create a harmonious transition to the diaphysis of the metatarsal. Capsule closure doubling the capsule to correct the position of the great toe and drawing the sesamoids medially. Skin suture. Dressing to splint the forefoot to support the corrected position of the great toe.

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Mobilisation with forefoot-relieving sandal for 6 weeks. If necessary, hallux valgus splint at night or adhesive strapping of the great toe for 6 weeks. Check radiograph before progressing to full loading of the forefoot and midfoot. The screws can be left in situ.

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Subluxations or dislocations of the MP joints II-V which can be reduced by a metatarsal shortening osteotomy. Excessive length of toes II-V e.g. after operation on the first ray. Contraindication: Destruction of the joint surface.

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
Skin incision on extensor aspect between two metatarsals commencing just distal to the MP joints (because of the intermetatarsal approach, the two adjacent joints can be reached from one incision). The extensor tendons are sought and held to one side, longitudinal opening of the joint capsule on the extensor aspect, plantarflexion of that toe, division of the lateral ligaments, on greater plantarflexion the metatarsal head emerges partially on the extensor aspect, pass two small Hohmann retractors around the exposed metatarsal head at the neck. Relatively shallow oblique osteotomy (roughly at the level of the foot) in the direction of the diaphysis commencing at the cartilage-bone junction (over a distance of approx. two centimetres). After this, the metatarsal head can be moved proximally over the desired distance. Removal of the bone edge projecting towards the joint on the extensor aspect, introduction of a self-tapping drill-fix-forefoot screw from the extensor aspect obliquely in the direction of the plantar joint surface. The tip of the screw must not emerge from the bone on the plantar aspect.

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Active exercises are possible immediately, if necessary adhesive strapping to maintain corrected position, forefoot-relieving shoe for 4 weeks. Check radiograph before progressing to full loading of the forefoot.
The screw can be left in situ.

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
Fracture compressing screw, diameter 4,0 mm / 3,0 mm, self-tapping, self-drilling

Bestell-Nr. Code N°	Länge Length	Anzahl Quantity
3.123.10	10 mm	3
3.123.12	12 mm	3
3.123.14	14 mm	3
3.123.16	16 mm	3
3.123.18	18 mm	3
3.123.20	20 mm	3
3.123.22	22 mm	3
3.123.24	24 mm	3
3.123.26	26 mm	3
3.123.28	28 mm	3
3.123.30	30 mm	3

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Drill-fix-forefoot screw, diameter 2,0 mm



Bestell-Nr. Code N°	Schraubenlänge screw length	Gesamtlänge total length	Anzahl Quantity
3.127.11	11 mm	26 mm	3
3.127.12	12 mm	27 mm	3
3.127.13	13 mm	28 mm	3
3.127.14	14 mm	29 mm	3

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
No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying, recording or otherwise without the prior written permission of the publisher.

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© Copyright Königsee Implantate GmbH
OSTEOSYNTHESE GmbH
ZUR OSTEOSYNTHESE GmbH**