

Unter Haßlen 20 · D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
e-mail: sales@rz-medizintechnik.com
http://www.rz-medizintechnik.com




FACELIFT

Atraumatic tissue
separation scissors



12/09



RZ Facelift Scissors Böttler

Appropriate instruments are essential for the separation of different tissue layers during facelift.

In particular an atraumatic separation according to the anatomical structures is highly important.

With many years of experience in the field of facelift surgery, Dr. Böttler has given the idea to create a new type of scissors: The front parts of the shanks have been reinforced by TC inserts. Most risks of skin perforation are remarkably reduced.

Major Advantages:

- successful and atraumatic separation of different tissue layers
- near tissues do not suffer any harm during separation
- less bleeding that enables a better healing



| 121-214-200

Faceliftscissors SC acc. to Böttler, reinforced, rounded shanks on the tip (acc. Kaye pattern), 20 cm



| 121-204-200

Faceliftscissors SC acc. to Böttler, reinforced, rounded shanks on the tip (acc. Sims pattern), 20 cm