



*progress in surgery*

Detachable  
NEEDLE  
HOLDERS




LAPAROSCOPY

# DETACHABLE DICK NEEDLE HOLDERS

This latest version of the RZ Laparoscopic Needle Holders convinces with advantages in cleaning purposes and with regard to a more cost-effective solution.


By detaching the needle holders in two or even three parts a higher grade of cleansing is achieved. The handle and outer shaft can be used each time while the insert can be changed. Even if inserts are worn-out, just the insert has to be replaced. Due to those aspects an economic solution in laparoscopic surgery is provided.




Detachable shaft from handle:  
convertible system, better  
cleaning possibilities



2-piece system:  
A extra thin Ø3mm System  
with 250mm and 330mm length  
for Pediatric and  
Mini-Laparoscopy




Secure screwing on the insert:  
Perfect fit of the insert and  
single cleaning  
(only 3-piece system)



Modular handles:  
Axial and straight handles  
available to fulfill the surgeon's  
needs and provide ergonomic  
comfort



3-piece system:  
Dick Needle holder can be  
fitted with different shafts and  
inserts, different lengths  
available



Anti-rotation feature:  
due to an extra rectangular  
connector, wrong insertion of  
the insert is excluded

In case of broken jaw, only the insert has to be replaced and the other parts of the instrument can continue to be used, which leads to a more cost-effective solution.

# DICK<sup>3</sup>

## 3-PIECE INSTRUMENT ∅ 5mm



Insert only	Complete instrument with axial handle	Complete instrument with straight handle	Complete instrument with angled handle
-------------	--	---	---



Straight jaw

330mm 300-405-130	300-480-630	300-480-130	300-480-160
450mm 300-406-130	300-480-730	300-480-230	300-480-260



Left curved jaw

330mm 300-405-140	300-480-640	300-480-140	300-480-170
450mm 300-406-140	300-480-740	300-480-240	300-480-270



Right curved jaw

330mm 300-405-141	300-480-641	300-480-141	300-480-171
450mm 300-406-141	300-480-741	300-480-241	300-480-271



Diamond dusted jaw, self-aligning

330mm 300-405-142	300-480-642	300-480-142	300-480-172
450mm 300-406-142	300-480-742	300-480-242	300-480-272






1x2 teeth jaw

330mm 300-405-131	300-480-631	300-480-131	300-480-161
450mm 300-406-131	300-480-731	300-480-231	300-480-261

## DICK<sup>2</sup>



### 2-PIECE INSTRUMENT ∅ 3mm

	Insert only	Complete instrument with straight handle	Complete instrument with angled handle
 Straight jaw	250mm 300-405-330 330mm 300-404-330	300-480-330 300-480-530	300-480-360 300-480-560
 Left curved jaw	330mm 300-405-340 450mm 300-404-340	300-480-340 300-480-540	300-480-370 300-480-570
 Right curved jaw	330mm 300-405-341 450mm 300-404-341	300-480-341 300-480-541	300-480-371 300-480-571

# SINGLE PARTS

## HANDLES

for combination with shaft and insert, two- and three-piece system



300-480-002

Axial Handle, strong, straight



300-480-000

Axial Handle, straight



300-480-001

Axial Handle, angled

## OUTER SHEATH

for combination with shaft and insert, three-pieces



300-480-033 330mm

300-480-045 450mm



*progress in surgery*

## RZ Medizintechnik GmbH

Unter Hasslen 20  
D-78532 Tuttlingen

Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50

[sales@rz-medizintechnik.com](mailto:sales@rz-medizintechnik.com)  
[www.rz-medizintechnik.com](http://www.rz-medizintechnik.com)

Media Code Number: 021490